
labscrip-utils

Release 3.1.0.dev80+g4581953

labscrip suite contributors

Nov 02, 2020

DOCUMENTATION

1 API Reference	3
1.1 Labscript Tools	3
1.2 Communications	9
1.3 GUI	11
1.4 Logging and Profiling	19
1.5 Module and File Tools	20
2 <i>labscript suite</i> components	25
Python Module Index	27
Index	29

labscrip-tutils contains code shared between multiple programs in the labscrip suite. This documentation is primarily for developers who might want to use this common code in their own custom features.

API REFERENCE

1.1 Labscript Tools

1.1.1 labscript_utils

```
labscript_utils.dedent(s)
```

Remove leading spaces from the first line of a string, all common leading indentation (spaces only) from subsequent lines, strip trailing spaces from all lines and replace single newlines prior to lines with the common indentation with spaces. Lines with additional indentation are kept verbatim. Good for unwrapping error messages etc that are in code as multiline triple-quoted strings.

```
labscript_utils.import_or_reload(modulename)
```

Behaves like ‘import modulename’ would, excepts forces the imported script to be rerun

1.1.2 labscript_utils.connections

```
class labscript_utils.connections.Connection(raw_row)
```

A class to represent a row in the connection table, present the contents as instance attributes after deserialising their contents, and providing default values for backward compatibility with older HDF5 files. Contains links to Connection objects for child devices of each device

```
_defaults = {'BLACS_connection': '', 'properties': {}, 'unit conversion class': None}
```

```
_deserialise(name, value)
```

deserialise one item of the row depending on what it is

```
_populate_relatives(table)
```

Populate child devices based on a list of other connection objects, and set self.parent to our parent device.

```
compare_to(other_connection)
```

```
diff(other)
```

```
find_by_name(name)
```

```
find_child(parent_name, parent_port)
```

```
print_details(indent)
```

```
property properties
```

```
property unit_conversion_params
```

```
class labscript_utils.connections.ConnectionTable(h5file, logging_prefix=None, exceptions_in_thread=False)
```

```
assert_superset (other)
compare_to (other)
find_by_name (name)
find_child (parent_name, parent_port)
getAttachedDevices ()
    Finds out which devices in the connection table are connected to BLACS, based on whether their 'BLACS_connection' attribute is non-empty. Returns a dictionary of them in the form {device_instance_name: labscript_class_name}

printDetails ()
removeDevice (device_name)
    Removes a device from the ConnectionTable, but keeps it in the raw_table. This can help make comparisons of connection tables fail for tables with broken devices.

labscript_utils.connections._ensure_str (s)
    convert bytestrings and numpy strings to python strings
```

1.1.3 labscript_utils.dict_diff

```
labscript_utils.dict_diff.dict_diff (dict1, dict2)
    Return the difference between two dictionaries as a dictionary of key: [val1, val2] pairs. Keys unique to either dictionary are included as key: [val1, '-'] or key: ['-', val2].
```

1.1.4 labscript_utils.labconfig

```
class labscript_utils.labconfig.EnvInterpolation
    Interpolation which expands environment variables in values, by post-filtering BasicInterpolation.before_get()

before_get (*args)

class labscript_utils.labconfig.LabConfig (config_path=PosixPath('/home/docs/labscript-
    suite/labconfig/build-12240137-project-
    608578-philipstarkey-labscript-utils.ini'),
    required_params=None, defaults=None)

exception NoOptionError (option, section)
    A requested option was not found.

exception NoSectionError (section)
    Raised when no section matches a requested option.

_abc_impl = <_abc_data object>

labscript_utils.labconfig.load_appconfig (filename)
    Load an .ini file and return a dictionary of its contents. All values will be converted to Python objects with ast.literal_eval(). All keys will be lowercase regardless of the written contents on the .ini file.

labscript_utils.labconfig.save_appconfig (filename, data)
    Save a dictionary as an ini file. The keys of the dictionary comprise the section names, and the values must themselves be dictionaries for the names and values within each section. All section values will be converted to strings with pprint.pformat().
```

1.1.5 labscript_utils.settings

```
class labscript_utils.settings.Settings (storage='hdf5',      file=None,      parent=None,
                                         page_classes=[])

add_settings_interface (setting_class)
close (*args, **kwargs)
create_dialog (goto_page=None)
get_value (settings_class, value_name)
load (name)
on_cancel (*args, **kwargs)
on_save (*args, **kwargs)
register_callback (callback)
remove_callback (callback)
```

1.1.6 labscript_utils.testing_utils

```
class labscript_utils.testing_utils.Any (types=<class 'object'>)
A class whose instances equal any object of the given type or tuple of types. For use with
mock.Mock.assert_called_with when you don't care what some of the arguments are

class labscript_utils.testing_utils.ThreadTestCase (*args, **kwargs)
Test case that runs tests in a new thread whilst providing a mainloop that allows running scripts in the current
thread. Those scripts can then be tested from the testing thread.

_mainloop ()
_run (*args, **kwargs)
    Called in a thread to run the tests
quit_mainloop ()
run (*args, **kwargs)
run_script_as_main (filepath)
static wait_for (condition_func, timeout=5, initial_poll_interval=0.005, max_poll_interval=0.5)
    Busy wait for a condition to be true. Uses exponential backoff so it's fast when things are fast and not a
    complete hog when they're not

class labscript_utils.testing_utils.dotdict
dot.notation access to dictionary attributes

class labscript_utils.testing_utils.monkeypatch (obj, name, mocked_attr)
Context manager to temporarily monkeypatch an object attribute with some mocked attribute
```

1.1.7 labscript_utils.properties

```
labscript_utils.properties._check_dicts(o)
labscript_utils.properties._decode_bytestrings(o)
    Decode all base64-encoded values (not keys) to bytestrings
labscript_utils.properties._default(o)
labscript_utils.properties._encode_bytestrings(o)
    Encode all bytestring values (not keys) to base64 with a prefix
labscript_utils.properties._get_con_table_properties(h5_file, device_name)
labscript_utils.properties._get_device_properties(h5_file, device_name)
labscript_utils.properties._get_unit_conversion_parameters(h5_file, device_name)
labscript_utils.properties.deserialise(value)
labscript_utils.properties.get(h5_file, device_name, location)
labscript_utils.properties.get_attribute(group, name)
    Return the attribute of the given name from the given HDF5 group, deserialising it if it has been encoded as JSON
labscript_utils.properties.get_attributes(group)
    Return attributes of a HDF5 group as a dict, deserialising any that have been encoded as JSON
labscript_utils.properties.is_json(value)
labscript_utils.properties.serialise(value)
labscript_utils.properties.set_attributes(group, attributes)
    Add attributes to a HDF5 group, serialising them to JSON if they do not map to native HDF5 datatypes
labscript_utils.properties.set_device_properties(h5_file, device_name, properties)
```

1.1.8 labscript_utils.unitconversions

```
labscript_utils.unitconversions.get_unit_conversion_class(fullname)
    import and return the unit conversion class with the given name. Ideally this is a fully qualified class name with an absolute import path, i.e. path.to.some.module.ClassName. But if it is just a single name, we fall back to looking through all classes defined in submodules. This allows backward compatibility with old shot files that do not have the full name saved.
```

Basic Unit Conversion Classes

```
class labscript_utils.unitconversions.test.test(calibration_parameters=None)
    Bases: labscript_utils.unitconversions.UnitConversionBase.UnitConversion

    A_from_base(volts)
    A_to_base(amps)
    Gauss_from_base(volts)
    Gauss_to_base(gauss)
    base_unit = 'MHz'
    derived_units = ['A', 'Gauss']
```

```

class labscript_utils.unitconversions.example.example1 (calibration_parameters=None)
Bases: labscript_utils.unitconversions.UnitConversionBase.UnitConversion

A_from_base (volts)
A_to_base (amps)
Gauss_from_base (volts)
Gauss_to_base (gauss)
base_unit = 'V'

class labscript_utils.unitconversions.example.example2 (calibration_parameters=None)
Bases: labscript_utils.unitconversions.UnitConversionBase.UnitConversion

base_unit = 'MHz'
detuned_MHz_from_base (mhz)
detuned_MHz_to_base (d_mhz)

class labscript_utils.unitconversions.example.example3 (calibration_parameters=None)
Bases: labscript_utils.unitconversions.UnitConversionBase.UnitConversion

W_from_base (vpp)
W_to_base (watts)
base_unit = 'Vpp'

class labscript_utils.unitconversions.UnitConversionBase.UnitConversion (params)
Bases: object

unit_list = {'G': 1000000000.0, 'M': 1000000.0, 'T': 1000000000000.0, 'k': 1000.0, 'm': 1000000000000000.0}
labscript_utils.unitconversions.UnitConversionBase.vectorise (method)

```

Unit Conversion Classes

```

class labscript_utils.unitconversions.aom.SineAom (calibration_parameters=None)
Bases: labscript_utils.unitconversions.NovaTechDDS9m.
NovaTechDDS9mAmpConversion

AOM calibration P(A) is very close to a sine for dipole trap AOM!

Power_from_base (amp)
Power_to_base (power)
base_unit = 'Arb'
fraction_from_base (amp)
fraction_to_base (fraction)

class labscript_utils.unitconversions.detuning.detuning (calibration_parameters=None)
Bases: labscript_utils.unitconversions.UnitConversionBase.UnitConversion

MHz_from_base (aom_frequency)
MHz_to_base (aom_frequency_MHz)
base_unit = 'Hz'
d_MHz_from_base (aom_frequency)

```

```
d_MHz_to_base (detuning_MHz)
derived_units = ['MHz', 'd_MHz', 'linewidths']

linewidths_from_base (aom_frequency)
linewidths_to_base (linewidths)

class labscript_utils.unitconversions.linear_coil_driver.BidirectionalCoilDriver (calibration_params)
Bases: labscript_utils.unitconversions.UnitConversionBase.UnitConversion

A_from_base (volts)
A_to_base (amps)
base_unit = 'V'
derived_units = ['A']

class labscript_utils.unitconversions.linear_coil_driver.UnidirectionalCoilDriver (calibration_params)
Bases: labscript_utils.unitconversions.linear_coil_driver.BidirectionalCoilDriver

A_from_base (volts)
A_to_base (amps)

class labscript_utils.unitconversions.NovaTechDDS9m.NovaTechDDS9mAmpConversion (calibration_params)
Bases: labscript_utils.unitconversions.UnitConversionBase.UnitConversion

base_unit = 'Arb'
hardware_from_base (arb)
hardware_to_base (hardware)

class labscript_utils.unitconversions.NovaTechDDS9m.NovaTechDDS9mFreqConversion (calibration_params)
Bases: labscript_utils.unitconversions.UnitConversionBase.UnitConversion

MHz_from_base (Hz)
MHz_to_base (MHz)
base_unit = 'Hz'

class labscript_utils.unitconversions.optotunelens.OptotuneLens (calibration_parameters=None)
Bases: labscript_utils.unitconversions.UnitConversionBase.UnitConversion

I_from_base (volts)
I_to_base (current)
base_unit = 'V'
derived_units = ['distance', 'I']
distance_from_base (volts)
distance_to_base (percentage)

class labscript_utils.unitconversions.quad_driver.quad_driver (calibration_parameters={'A_min':
- 0.09, 'A_offset':
- 0.642724,
'A_per_V':
19.9757,
'Gcm_per_A':
1.88679})
Bases: labscript_utils.unitconversions.UnitConversionBase.UnitConversion
```

```

A_from_base(volts)
A_to_base(arg)
Gcm_from_base(volts)
Gcm_to_base(gauss_per_cm)
base_unit = 'V'
derived_units = ['A', 'Gcm']

class labscript_utils.unitconversions.quad_monitor.quad_monitor(calibration_parameters={'A_offset': 0.043200000000000016, 'A_per_V': 20.032, 'Gcm_per_A': 1.88679})

Bases: labscript_utils.unitconversions.UnitConversionBase.UnitConversion

A_from_base(volts)
A_to_base(amps)
Gcm_from_base(volts)
Gcm_to_base(gauss_per_cm)
base_unit = 'V'
derived_units = ['A', 'Gcm']

```

1.2 Communications

1.2.1 labscript_utils.h5_lock

```

class labscript_utils.h5_lock.File(name, mode=None, driver=None, libver=None, **kwds)

    _abc_impl = <_abc_data object>
    close()
        Close the file. All open objects become invalid
    labscript_utils.h5_lock.hack_locks_onto_h5py()

```

1.2.2 labscript_utils.ls_zprocess

```

class labscript_utils.ls_zprocess.Context(io_threads=1, shared_secret=None)
    Subclass of zprocess.security.SecureContext configured with settings from labconfig, substitutable for a zmq.Context. Can be instantiated to get a unique context, or call the .instance() classmethod to possibly get an already-existing one. Only use the latter if you do not intend to terminate the context.

    classmethod instance()
        Returns a shared instance with the same shared secret, if there is one, otherwise creates it. If an instance already exists, io_threads will be ignored, otherwise it will be used in the new instance. Takes into account subclasses such that a subclass calling this method will always get back an instance of its own class

```

socket(*args, **kwargs)

Create a Socket associated with this Context.

Parameters

- **socket_type** (*int*) – The socket type, which can be any of the 0MQ socket types: REQ, REP, PUB, SUB, PAIR, DEALER, ROUTER, PULL, PUSH, etc.
- **kwargs** – will be passed to the `__init__` method of the socket class.

`labsclient_utils.ls_zprocess.Event(*args, **kwargs)``labsclient_utils.ls_zprocess.Handler(*args, **kwargs)``labsclient_utils.ls_zprocess.Lock(*args, **kwargs)``class labsclient_utils.ls_zprocess.ProcessTree(shared_secret=None, al-
low_insecure=False, zlock_host=None,
zlock_port=7339, zlog_host=None,
zlog_port=7340)`

A singleton zprocess.ProcessTree configured with settings from labconfig for security, zlock and zlog. Being a singleton is not enforced - the class can still be instantiated as normal - but calling the `.instance()` classmethod will give the singleton.

```
_instance = <labsclient_utils.ls_zprocess.ProcessTree object>  
  
classmethod instance()  
  
labsclient_utils.ls_zprocess.RemoteProcessClient(host, port=None)  
  
class labsclient_utils.ls_zprocess.ZMQClient  
A singleton zprocess.ZMQClient configured with settings from labconfig for security. Being a singleton is not enforced - the class can still be instantiated as normal - but calling the .instance() classmethod will give the singleton.  
  
_instance = None  
  
classmethod instance()  
  
class labsclient_utils.ls_zprocess.ZMQServer(port=None, dtype='pyobj', pull_only=False,  
bind_address='tcp://0.0.0.0', time-  
out_interval=None, **kwargs)
```

A ZMQServer configured with security settings from labconfig

`labsclient_utils.ls_zprocess.connect_to_zlock_server()``labsclient_utils.ls_zprocess.ensure_connected_to_zlog()`

Ensure we are connected to a zlog server. If one is not running and we are the top-level process, start a zlog server configured according to LabConfig.

`labsclient_utils.ls_zprocess.get_config()`

Get relevant options from LabConfig, substituting defaults where appropriate and return as a dict

`labsclient_utils.ls_zprocess.zmq_get(*args, **kwargs)``labsclient_utils.ls_zprocess.zmq_get_multipart(*args, **kwargs)``labsclient_utils.ls_zprocess.zmq_get_raw(*args, **kwargs)``labsclient_utils.ls_zprocess.zmq_get_string(*args, **kwargs)``labsclient_utils.ls_zprocess.zmq_push(*args, **kwargs)``labsclient_utils.ls_zprocess.zmq_push_multipart(*args, **kwargs)``labsclient_utils.ls_zprocess.zmq_push_raw(*args, **kwargs)`

```
labscript_utils.ls_zprocess.zmq_push_string(*args, **kwargs)
```

1.2.3 labscript_utils.remote

Script to run a zprocess.remote server configured according to LabConfig. Run with:

```
python -m labscript_utils.remote [--daemon] [--no-tui]
```

If --daemon is specified, the server will be started in the background. If --no-tui is specified, the server will run with ordinary terminal output instead of with the interactive text-based user interface (TUI).

```
labscript_utils.remote.main()
```

1.2.4 labscript_utils.shared_drive

```
labscript_utils.shared_drive.path_to_agnostic(path)
labscript_utils.shared_drive.path_to_local(path)
```

1.2.5 labscript_utils.zlock

Script to run a zlock server configured according to LabConfig. Run with:

```
python -m labscript_utils.zlock [--daemon]
```

If --daemon is specified, the zlock server will be started in the background.

```
labscript_utils.zlock.main()
```

1.2.6 labscript_utils.zlog

Script to run a zlog server configured according to LabConfig. Run with:

```
python -m labscript_utils.zlog [--daemon]
```

If --daemon is specified, the zlog server will be started in the background.

```
labscript_utils.zlog.main()
```

1.3 GUI

1.3.1 labscript_utils.splash

```
class labscript_utils.splash.Splash(imagepath)

    BG = '#ffffff'
    alpha = 0.875
    h = 230
    icon_frac = 0.65
```

```
imheight = 150
imwidth = 150
paintEvent (self, QPaintEvent)
show (self)
update_text (text)
w = 250
```

1.3.2 labscrip_utils.qtwidgets

Helper Widgets

```
class labscrip_utils.qtwidgets.dragdroptab.DragDropTabBar (parent, group_id)
Bases: labscrip_utils.qtwidgets.dragdroptab._BaseDragDropTabBar

FLUSH_GAP = 5
SCROLL_BUTTON_GAP = 2
SCROLL_BUTTON_WIDTH = 15
property drag_in_progress
property dragged_tab_index
property dragged_tab_parent
ensure_visible (index, prefer_left=True)
insertion_index_at (pos)
    Compute at which index the tab with given upper left corner position in global coordinates should be inserted into the tabBar.
is_dragged_tab (index)
    Return whether the tab at the given index is currently being dragged
limbo = None
minimumSizeHint (self) → QSize
mouseMoveEvent (event)
    Update the parent of the tab to be the DragDropTabWidget under the mouse, if any, otherwise update it to the limbo object. Update the position of the tab in the widget it's in.
mousePressEvent (event)
    Take note of the tab that was clicked so it can be dragged on mouseMoveEvents
mouseReleaseEvent (event)
    Same as mouseMove event - update the DragDropTabWidget and position of the tab to the current mouse position. Unless the mouse position is outside of any widgets at the time of mouse release, in which case move the tab to its last known parent and position.
on_scroll_button_clicked (button)
paintEvent (self, QPaintEvent)
paint_tab (index, painter, option)
setUsesScrollButtons (self, bool)
```

set_tab_parent (*dest, index=None, pos=None*)
 Move the tab to the given parent DragDropTabBar if it's not already there. *index=None* will determine the insertion index from the given mouse position.

sizeHint (*self*) → **QSize**

tabAt (*self, QPoint*) → **int**

tabInserted (*self, int*)

tabLayoutChange (*self*)

tabRect (*self, int*) → **QRect**

tabRemoved (*self, int*)

tab_widgets = {}

update (*self*)

 update(*self, QRect*) update(*self, QRegion*) update(*self, int, int, int, int*)

update_dragged_tab_animation_pos (*pos*)

update_scroll_button_state ()

update_tab_index (*index, pos*)

 Check if the tab at the given index, being dragged by the mouse at the given position, needs to be moved.
 Move it and return the new index.

widgetAt (*pos*)

 If the given position is over a DragDropTabBar belonging to the current group, return the DragDropTabBar.
 If it is over a TabWidget in the same group that has no tabs, or the dragged tab as its only tab, return its DragDropTabBar. Otherwise return the limbo object.

class labscript_utils.qtwidgets.dragdroptab.**DragDropTabWidget** (*group_id=None, ac-cept_drops_bar_only=False*)

Bases: **PyQt5.QtWidgets.QTabWidget**

A tab widget that supports dragging and dropping of tabs between tab widgets that share a group_id. A group_id of None indicates that tab dragging is disabled.

setElideMode (*self, Qt.TextElideMode*)

setUsesScrollButtons (*self, bool*)

class labscript_utils.qtwidgets.dragdroptab.**Tab** (*widget, text, data, text_color, tooltip, whats_this, button_left, button_right, icon*)

Bases: **tuple**

property button_left

 Alias for field number 6

property button_right

 Alias for field number 7

property data

 Alias for field number 2

property icon

 Alias for field number 8

property text

 Alias for field number 1

```
property text_color
    Alias for field number 3

property tooltip
    Alias for field number 4

property whats_this
    Alias for field number 5

property widget
    Alias for field number 0

class labscrip_utils.qtwidgets.dragdroptab.TabAnimation(parent)
Bases: PyQt5.QtCore.QAbstractAnimation

animate_limbo(limbo, index)
    If the floating tab in limbo is being sucked back into one of our tabs, then we can animate that by hiding
    the relevant tab rect off to the side somewhere whilst the floating tab swoops in.

duration(self) → int
ensure_running()
on_tab_moved(source_index, dest_index)
tabInserted(index)
tabRemoved(index)
target(i)
    Return the target position we are animating toward for a tab
tau = 60.0
updateCurrentTime(self, int)

class labscrip_utils.qtwidgets.dragdroptab.debug
Bases: object

DEBUG = False
depth = 0

classmethod trace(f)
    decorator to print function entries and exits

class labscrip_utils.qtwidgets.elide_label.ElideScrollArea(*args, **kwargs)
Bases: PyQt5.QtWidgets.QScrollArea

A ScrollArea for containing a label that we want to elide. The elision is attained by just letting the text we don't
want to see be scrolled off to the side with the scrollbars hidden.

event(self, QEvent) → bool
minimumSizeHint(self) → QSize
setElideMode(elideMode)
setWidget(self, QWidget)
sizeHint(self) → QSize

class labscrip_utils.qtwidgets.elide_label.ElidedLabelContainer(label)
Bases: PyQt5.QtWidgets.QWidget
```

A QWidget to contain a QLabel with a single line of (possibly rich) text that we want to elide. The elision is obtained by putting the QLabel in a QScrollArea and having the QScrollarea only show the part of the text we want to see. An extra label with the elision indication “...” is also inserted next to the QScrollArea.

```
elideMode ()  

event (self, QEvent) → bool  

minimumSizeHint (self) → QSize  

setElideMode (elideMode)  

sizeHint (self) → QSize  

update_elide_widget ()
```

`labscrip_utils.qtwidgets.elide_label.elide_label(label, layout, elide_mode)`

Take an existing label that is in a layout, and wrap it in our widgets that elide the text, and insert it back into the layout. This is a hack that allows us to elide a QLabel with a single line of (possibly rich) text, a task that seems pretty much impossible to do in any kosher way.

This function is for modifying an existing label already in a layout, but if you are programmatically creating a label, then you can wrap it in ElidedLabelContainer(label) before inserting it into a layout or other container widget, which is more flexible than this function which only works if the label is in a QHBoxLayout

```
class labscrip_utils.qtwidgets.fingertab.FingerTabBarWidget (parent=None,  
 *args, **kwargs)
```

Bases: PyQt5.QtWidgets.QTabBar

paintEvent (*self*, *QPaintEvent*)

tabSizeHint (*self*, *int*) → *QSize*

```
class labscrip_utils.qtwidgets.fingertab.FingerTabWidget (parent, *args)
```

Bases: PyQt5.QtWidgets.QTabWidget

A QTabWidget equivalent which uses our FingerTabBarWidget

```
class labscrip_utils.qtwidgets.headerview_with_widgets.HorizontalHeaderViewWithWidgets (mo  
 par  
 ent:
```

Bases: PyQt5.QtWidgets.QHeaderView

A QHeaderView that supports inserting arbitrary widgets into sections. Use setWidget(logical_index, widget) to set and setWidget(logical_index, None) to unset. Decorations, checkboxes or anything other than text in the headers containing widgets is unsupported, and may result in garbled output

do_update_widget_positions ()

eventFilter (*target*, *event*)

Ensure we don't leave the cursor set as a resize handle when the mouse moves onto a child widget:

hideSection (*self*, *int*)

on_columnsInserted (*parent*, *logical_first*, *logical_last*)

on_columnsRemoved (*parent*, *logical_first*, *logical_last*)

sectionSizeFromContents (*self*, *int*) → *QSize*

setSectionHidden (*self*, *int*, *bool*)

setStyleSheet (*self*, *str*)

setWidget (*logical_index*, *widget=None*)

showEvent (*self*, *QShowEvent*)

```
showSection (self, int)
stylesheet = '\n QHeaderView::section {\n /* Will be set dynamically: */\n padding-top: 0px;\n padding-bottom: 0px;\n border-bottom: 1px solid black;\n }'
thinspace = '\u2009'
update_indent()
update_widget_positions()
viewportEvent (self, QEvent) → bool

class labscript_utils.qtwidgets.InputPlotWindow.PlotWindow (*args, **kwargs)
Bases: zprocess.process_tree.Process

run (connection_name, hardware_name, device_name)
    The method that gets called in the subprocess. To be overridden by subclasses

setTopLevelWindow()

update_plot (new_data)

class labscript_utils.qtwidgets.outputbox.OutputBox (container, scroll-  
back_lines=1000)
Bases: qtutils.outputbox.OutputBox

A subclass of qtutils.outputbox.OutputBox configured with security from labconfig.

class labscript_utils.qtwidgets.toolpalette.ToolPalette (parent, name, *args,  
**kwargs)
Bases: PyQt5.QtWidgets.QScrollArea

addWidget (widget, force_relayout=True)
insertWidget (index, widget, force_relayout=True)
minimumSize (self) → QSize
minimumSizeHint (self) → QSize
resizeEvent (self, QResizeEvent)
sizeHint (self) → QSize
updateMinimumSize()

class labscript_utils.qtwidgets.toolpalette.ToolPaletteGroup (*args, **kwargs)
Bases: PyQt5.QtWidgets.QVBoxLayout

add_to_linked_width_group (width_group_name, name)
append_new_palette (name, *args, **kwargs)
create_linked_width_group (width_group_name, names)
get_index_from_name (name)
get_name_from_index (index)
get_palette (name)
get_palette_by_index (index)
has_palette (name)
hide_palette (name)
hide_palette_by_index (index)
insert_new_palette (index, name, *args, **kwargs)
```

```

remove (name)
remove_by_index (index)
remove_from_linked_width_group (width_group_name, name)
reorder_palette (name, new_index)
reorder_palette_by_index (old_index, new_index)
show_palette (name)
show_palette_by_index (index)
property widths_linked

```

Input/Output Widgets

```

class labscript_utils.qtwidgets.analoginput.AnalogInput (device_name,
                                                       hardware_name,
                                                       connection_name='-',
                                                       display_name=None, horizontal_alignment=False,
                                                       parent=None)

Bases: PyQt5.QtWidgets.QWidget

get_AI ()
open_plot_window ()
set_AI (AI, notify_old_AI=True, notify_new_AI=True)
set_value (value)

class labscript_utils.qtwidgets.analogoutput.AnalogOutput (hardware_name,
                                                               connection_name='-',
                                                               display_name=None,
                                                               horizontal_alignment=False,
                                                               parent=None)

Bases: PyQt5.QtWidgets.QWidget

block_combobox_signals ()
block_spinbox_signals ()
connect_value_change (func)
disconnect_value_change ()
eventFilter (self, QObject, QEvent) → bool
get_AO ()
lock (notify_ao=True)
property selected_unit
set_AO (AO, notify_old_AO=True, notify_new_AO=True)
set_combobox_model (model)
set_limits (lower, upper)
set_num_decimals (decimals)

```

```
    set_selected_unit (unit)
    set_spinbox_value (value, unit)
    set_step_size (step)
    unblock_combobox_signals ()
    unblock_spinbox_signals ()
    unlock (notify_ao=True)

class labscript_utils.qtwidgets.analogoutput.NoStealFocusDoubleSpinBox (*args,
                                                               **kwargs)
Bases: PyQt5.QtWidgets.QDoubleSpinBox
A QDoubleSpinBox that doesn't steal focus as you scroll over it with a mouse wheel.

    focusInEvent (self, QFocusEvent)
    focusOutEvent (self, QFocusEvent)
    wheelEvent (self, QWheelEvent)

class labscript_utils.qtwidgets.ddsoutput.DDSSOutput (hardware_name,
                                                       connection_name='-',      parent=None)
Bases: PyQt5.QtWidgets.QWidget
    get_sub_widget (subchnl)
    hide_sub_widget (subchnl)
    show_sub_widget (subchnl)

class labscript_utils.qtwidgets.digitaloutput.DigitalOutput (*args, **kwargs)
Bases: PyQt5.QtWidgets.QPushButton
    eventFilter (self, QObject, QEvent) → bool
    get_DO ()
    lock (notify_do=True)
    set_DO (DO, notify_old_DO=True, notify_new_DO=True)
    property state
    unlock (notify_do=True)

class labscript_utils.qtwidgets.digitaloutput.InvertedDigitalOutput (*args,
                                                               **kwargs)
Bases: labscript_utils.qtwidgets.digitaloutput.DigitalOutput
    property state

class labscript_utils.qtwidgets.enumoutput.EnumOutput (hardware_name,
                                                       connection_name='-',      display_name=None,      horizontal_alignment=False,
                                                       parent=None)
Bases: PyQt5.QtWidgets.QWidget
    block_combobox_signals ()
    connect_value_change (func)
    disconnect_value_change ()
```

```

eventFilter(self, QObject, QEvent) → bool
get_EO()
lock(notify_eo=True)
property selected_index
property selected_option
set_EO(EO, notify_old_EO=True, notify_new_EO=True)
set_combobox_model(model)
unblock_combobox_signals()
unlock(notify_eo=True)

class labscript_utils.qtwidgets.imageoutput.BrowseButton(image_output,      *args,
                                                               **kwargs)
Bases: PyQt5.QtWidgets.QPushButton

browse()

eventFilter(self, QObject, QEvent) → bool

class labscript_utils.qtwidgets.imageoutput.ImageOutput(name, width, height, *args,
                                                               **kwargs)
Bases: PyQt5.QtWidgets.QWidget

eventFilter(self, QObject, QEvent) → bool

get_Image()

imageUpdated

lock(notify_Image=True)

set_Image(Image, notify_old_Image=True, notify_new_Image=True)
unlock(notify_Image=True)

property value

class labscript_utils.qtwidgets.imageoutput.ImageView(*args, **kwargs)
Bases: PyQt5.QtWidgets.QGraphicsView

contextMenuEvent(self, QContextMenuEvent)

```

1.4 Logging and Profiling

1.4.1 labscript_utils.impprof

```

class labscript_utils.impprof._ProfilingImporter

disable()
enable(threshold=0.1)
profiling_import(name, *args, **kwargs)

labscript_utils.impprof.disable()
labscript_utils.impprof.enable(threshold=0.1)

```

1.4.2 labscrip_utils.memprof

```
class labscrip_utils.memprof.MemoryProfiler
    Class to count number instances of each type in the interpreter in order to detect Python memory leaks

    check()
    count_types()
    start(filepath='memprof.txt')
    write_to_file(types)

labscrip_utils.memprof.check()
labscrip_utils.memprof.start(filepath='memprof.txt')
```

1.4.3 labscrip_utils.setup_logging

```
class labscrip_utils.setup_logging.LessThanFilter(less_than)

filter(record)
    Determine if the specified record is to be logged.

    Is the specified record to be logged? Returns 0 for no, nonzero for yes. If deemed appropriate, the record
    may be modified in-place.

labscrip_utils.setup_logging.setup_logging(program_name, log_level=10, terminal_level=20, maxBytes=52428800, backupCount=1)
```

1.4.4 labscrip_utils.tracelog

```
labscrip_utils.tracelog.log(log_path=None, module_names=, sub=False, all=False,
```

Trace and log Python execution.

output includes the time, thread name, containing function name, line number and source line. Indentation before the thread name represents stack depth, indentation before source line is as in the source line itself.

log_path: the path of the desired output file to write to, or None for stdout (default=None) module_names: list of module names that tracing is desired for (default=()) sub: whether submodules of the above modules should be traced (default=False) all: whether all modules should be traced, in which case module_names is ignored (default=False) mode: mode to open the output file in, if log_path is not None (default='w')

1.5 Module and File Tools

1.5.1 labscrip_utils.double_import_denier

```
class labscrip_utils.double_import_denier.DoubleImportDenier
```

A module finder that tracks what's been imported and disallows multiple imports of the same module under different names, raising an exception upon detecting that this has occurred

```
_format_tb(tb)
```

Take a formatted traceback as returned by traceback.format_stack() and remove lines that are solely about us and the Python machinery, leaving only lines pertaining to the user's code

```

    _raise_error(path, name, tb, other_name, other_tb)
    _restore_tracebacklimit_after_exception()
        Record the current value of sys.tracebacklimit, if any, and set a temporary sys.excepthook to restore it to
        that value (or delete it) after the next exception.

    find_spec(fullname, path=None, target=None)

labscript_utils.double_import_denier.disable()
labscript_utils.double_import_denier.enable()

```

1.5.2 labscript_utils.filewatcher

```

class labscript_utils.filewatcher.FileWatcher(callback, files=None, folders=None,
                                              clean_modified_info=None,          hash-
                                              able_types=None, interval=1, **kwargs)

    _modified_info_of_file(name)
    add_file(path)
    add_files(files, clean_modified_info=None)
    add_folder(folder)
    add_folders(folders, clean_modified_info=None)
    check(trigger_callback=True)
    get_clean_modified_info()
    get_modified_info()
    get_modified_times()
    mainloop()
    stop()
    update_files(folders=None, trigger_callback=True, recursive=True)
        Refresh the watchlist of files (FileWatcher.files) by checking the folders kwarg or Filewatcher.folders if
        this is not specified.

```

1.5.3 labscript_utils.modulewatcher

```

class labscript_utils.modulewatcher.ModuleWatcher(debug=False)
    A watcher that reloads modules that have been modified on disk

    Only reloads modules imported after instantiation. Does not reload C extensions.

    Parameters debug (bool, optional) – When True, prints debugging information when
        reloading modules.

    check()
    mainloop()
    unload()

```

1.5.4 labscrip_utils.versions

```
exception labscrip_utils.versions.BrokenInstall
```

```
class labscrip_utils.versions.NoVersionInfo
```

```
class labscrip_utils.versions.NotFound
```

```
exception labscrip_utils.versions.VersionException
```

```
labscrip_utils.versions._get_literal_version(filename)
```

Tokenize a source file and return any __version__ = <version> literal defined in it.

Parameters `filename` (`str`) – The path to the file to tokenize.

Returns Any version literal found matching the above criteria, or None.

```
labscrip_utils.versions._get_metadata_version(project_name, import_path)
```

Gets the package metadata version.

Parameters

- `project_name` (`str`) – The package name (e.g. the name used when pip installing the package).

- `import_path` (`str`) – The path to the folder containing the installed package.

Raises `BrokenInstall` – Raised if the package installation is corrupted (multiple packages matching the given arguments were found). May occur if (un)installation for a particular package version only partially completed.

Returns The metadata version for a package with the given project name located at the given import path, or None if there is no such package.

```
labscrip_utils.versions.check_version(module_name, at_least, less_than, version=None, project_name=None)
```

Checks if a module version is within specified bounds.

Checks that the version of the given module is at least and less than the given version strings. This function uses `get_version()` to determine version numbers without importing modules. In order to do this, `project_name` must be provided if it differs from `module_name`. For example, pyserial is imported as ‘serial’, but the project name, as passed to a ‘pip install’ command, is ‘pyserial’. Therefore to check the version of pyserial, pass in `module_name='serial'` and `project_name='pyserial'`. You can also pass in a version string yourself, in which case no inspection of packages will take place.

Parameters

- `module_name` (`str`) – The name of the module to check.
- `at_least` (`str`) – The minimum acceptable module version.
- `less_than` (`str`) – The minimum unacceptable module version. Usually this would be the next major version if the package follows semver.
- `version` (`str`, optional) – The current version of the installed package. Useful when the package version is stored in a non-standard location.
- `project_name` (`str`, optional) – The package name (e.g. the name used when pip installing the package). This must be specified if it does not match the module name.

Raises `VersionException` – if the module was not found or its version could not be determined.

```
labscrip_utils.versions.get_import_path(import_name)
```

Get which entry in sys.path a module would be imported from, without importing it.

Parameters `import_name` (`str`) – The module name.

Raises

- **ModuleNotFoundError** – Raised if the module is not installed.
- **NotImplementedError** – Raised if the module is a “namespace package”. Support for namespace packages is not currently available.

Returns The path to the folder containing the module.

Return type str

```
labscript_utils.versions.get_version(import_name,          project_name=None,      im-
                                     port_path=None)
```

Try very hard to get the version of a package without importing it.

If import_path is not given, first find where it would be imported from, without importing it. Then look for metadata in the same import path with the given project name (note: this is not always the same as the import name, it is the name for example you would ask pip to install). If that is found, return the version info from it. Otherwise look for a `__version__.py` file in the package directory, or a `__version__ = <version>` literal defined in the package source (without executing it).

Parameters

- **import_name** (str) – The module name.
- **project_name** (str, optional) – The package name (e.g. the name used when pip installing the package). This must be specified if it does not match the module name.
- **import_path** (str, optional) – The path to the folder containing the installed package.

Raises **NotImplementedError** – Raised if the module name contains a period. Only top-level packages are supported at this time.

Returns The version literal of the package. If the package cannot be found, `NotFound` is returned.

If the version cannot be obtained in the above way, or if the version was found but was None, `NoVersionInfo` is returned.

CHAPTER
TWO

LABSCRIPT SUITE COMPONENTS

The *labscript suite* is modular by design, and is comprised of:

Table 1: Python libraries

	labscript — Expressive composition of hardware-timed experiments
	labscript-devices — Plugin architecture for controlling experiment hardware
	labscript-utils — Shared modules used by the <i>labscript suite</i>

Table 2: Graphical applications

	runmanager — Graphical and remote interface to parameterized experiments
	blacs — Graphical interface to scientific instruments and experiment supervision
	lyse — Online analysis of live experiment data
	runviewer — Visualize hardware-timed experiment instructions

PYTHON MODULE INDEX

|

labscript_utils, 3
labscript_utils.connections, 3
labscript_utils.dict_diff, 4
labscript_utils.double_import_denier,
 20
labscript_utils.filewatcher, 21
labscript_utils.h5_lock, 9
labscript_utils.impprof, 19
labscript_utils.labconfig, 4
labscript_utils.ls_zprocess, 9
labscript_utils.memprof, 20
labscript_utils.modulewatcher, 21
labscript_utils.properties, 6
labscript_utils.qtwidgets, 12
labscript_utils.qtwidgets.analoginput,
 17
labscript_utils.qtwidgets.analogoutput,
 17
labscript_utils.qtwidgets.ddsoutput, 18
labscript_utils.qtwidgets.digitaloutput,
 18
labscript_utils.qtwidgets.dragdroptab,
 12
labscript_utils.qtwidgets.elide_label,
 14
labscript_utils.qtwidgets.enumoutput,
 18
labscript_utils.qtwidgets.fingertab, 15
labscript_utils.qtwidgets.headerview_with_widgets,
 15
labscript_utils.qtwidgets.imageoutput,
 19
labscript_utils.qtwidgets.InputPlotWindow,
 16
labscript_utils.qtwidgets.outputbox, 16
labscript_utils.qtwidgets.palette,
 16
labscript_utils.remote, 11
labscript_utils.settings, 5
labscript_utils.setup_logging, 20
labscript_utils.shared_drive, 11
labscript_utils.splash, 11
labscript_utils.testing_utils, 5
labscript_utils.tracelog, 20
labscript_utils.unitconversions, 6
labscript_utils.unitconversions.aom, 7
labscript_utils.unitconversions.detuning,
 7
labscript_utils.unitconversions.example,
 6
labscript_utils.unitconversions.linear_coil_driver,
 8
labscript_utils.unitconversions.NovaTechDDS9m,
 8
labscript_utils.unitconversions.optotunelens,
 8
labscript_utils.unitconversions.quad_driver,
 8
labscript_utils.unitconversions.quad_monitor,
 9
labscript_utils.unitconversions.test, 6
labscript_utils.unitconversions.UnitConversionBase,
 7
labscript_utils.versions, 22
labscript_utils.zlock, 11
labscript_utils.zlog, 11

INDEX

Symbols

_ProfilingImporter (class in labscript_utils.impprof), 19
_abc_impl (labscript_utils.h5_lock.File attribute), 9
_abc_impl (labscript_utils.labconfig.LabConfig attribute), 4
_check_dicts () (in module labscript_utils.properties), 6
_decode_bytestrings () (in module labscript_utils.properties), 6
_default () (in module labscript_utils.properties), 6
_defaults (labscript_utils.connections.Connection attribute), 3
_deserialise () (labscript_utils.connections.Connection method), 3
_encode_bytestrings () (in module labscript_utils.properties), 6
_ensure_str () (in module labscript_utils.connections), 4
_format_tb () (labscript_utils.double_import_denier.DoubleImportDenier method), 20
_get_con_table_properties () (in module labscript_utils.properties), 6
_get_device_properties () (in module labscript_utils.properties), 6
_get_literal_version () (in module labscript_utils.versions), 22
_get_metadata_version () (in module labscript_utils.versions), 22
_get_unit_conversion_parameters () (in module labscript_utils.properties), 6
_instance (labscript_utils.ls_zprocess.ProcessTree attribute), 10
_instance (labscript_utils.ls_zprocess.ZMQClient attribute), 10
_mainloop () (labscript_utils.testing_utils.ThreadTestCase method), 5
_modified_info_of_file () (labscript_utils.filewatcher.FileWatcher method), 21
_populate_relatives () (labscript_utils.connections.Connection method), 3
_raise_error () (labscript_utils.double_import_denier.DoubleImportDenier method), 21
_restore_tracebacklimit_after_exception () (labscript_utils.double_import_denier.DoubleImportDenier method), 21
_run () (labscript_utils.testing_utils.ThreadTestCase method), 5

A

A_from_base () (labscript_utils.unitconversions.example.example1 method), 7
A_from_base () (labscript_utils.unitconversions.linear_coil_driver.BidirectionalCoilDriver method), 8
A_from_base () (labscript_utils.unitconversions.linear_coil_driver.UnidirectionalCoilDriver method), 8
A_from_base () (labscript_utils.unitconversions.quad_driver.quad_driver method), 8
A_from_base () (labscript_utils.unitconversions.quad_monitor.quad_monitor method), 9
A_from_base () (labscript_utils.unitconversions.test.test method), 6
A_to_base () (labscript_utils.unitconversions.example.example1 method), 7
A_to_base () (labscript_utils.unitconversions.linear_coil_driver.BidirectionalCoilDriver method), 8
A_to_base () (labscript_utils.unitconversions.linear_coil_driver.UnidirectionalCoilDriver method), 8
A_to_base () (labscript_utils.unitconversions.quad_driver.quad_driver method), 9
A_to_base () (labscript_utils.unitconversions.quad_monitor.quad_monitor method), 9
A_to_base () (labscript_utils.unitconversions.test.test method), 6

```

        method), 6
add_file() (labscript_utils.filewatcher.FileWatcher
            method), 21
add_files() (labscript_utils.filewatcher.FileWatcher
            method), 21
add_folder() (lab-
            script_utils.filewatcher.FileWatcher    method), 21
add_folders() (lab-
            script_utils.filewatcher.FileWatcher    method), 21
add_settings_interface() (lab-
            script_utils.settings.Settings method), 5
add_to_linked_width_group() (lab-
            script_utils.qtwidgets.toolpalette.ToolPaletteGroup
            method), 16
addWidget () (labscript_utils.qtwidgets.toolpalette.ToolPalette
            method), 16
alpha (labscript_utils.splash.Splash attribute), 11
AnalogInput (class in lab-
            script_utils.qtwidgets.analoginput), 17
AnalogOutput (class in lab-
            script_utils.qtwidgets.analogoutput), 17
animate_limbo() (lab-
            script_utils.qtwidgets.dragdroptab.TabAnimation
            method), 14
Any (class in labscript_utils.testing_utils), 5
append_new_palette() (lab-
            script_utils.qtwidgets.toolpalette.ToolPaletteGroup
            method), 16
assert_superset () (lab-
            script_utils.connections.ConnectionTable
            method), 3
base_unit (labscript_utils.unitconversions.aom.SineAom
            check () (in module labscript_utils.memprof), 20
            attribute), 7
base_unit (labscript_utils.unitconversions.detuning.detuning
            method), 21
            attribute), 7
base_unit (labscript_utils.unitconversions.example.example1
            method), 20
            attribute), 7
base_unit (labscript_utils.unitconversions.example.example2
            method), 21
            attribute), 7
base_unit (labscript_utils.unitconversions.example.example3
            method), 21
            attribute), 7
base_unit (labscript_utils.unitconversions.linear_coil_driver.BidirectionalCoilDriver
            BiDirectionalCoilDriverSettings.Settings method), 5
            attribute), 8
base_unit (labscript_utils.unitconversions.NovaTechDDS9m.NovaTechDDS9m
            NovaTechDDS9mAmpConvex.Connection method),
            attribute), 8
base_unit (labscript_utils.unitconversions.NovaTechDDS9m.NovaTechDDS9m
            NovaTechDDS9mFreqConversion
            attribute), 8
base_unit (labscript_utils.unitconversions.optotunelens.OptotuneLens
            method), 4
            attribute), 8
base_unit (labscript_utils.unitconversions.quad_driver.quad_driver
            base_unit (labscript_utils.unitconversions.quad_monitor.quad_monitor
            attribute), 9
            attribute), 9
base_unit (labscript_utils.unitconversions.test.test_at-
            tribute), 6
before_get () (lab-
            script_utils.labconfig.EnvInterpolation
            method), 4
BG (labscript_utils.splash.Splash attribute), 11
BidirectionalCoilDriver (class in lab-
            script_utils.unitconversions.linear_coil_driver),
            8
block_combobox_signals() (lab-
            script_utils.qtwidgets.analogoutput.AnalogOutput
            method), 17
k_combobox_signals() (lab-
            script_utils.qtwidgets.enumoutput.EnumOutput
            method), 18
block_spinbox_signals() (lab-
            script_utils.qtwidgets.analogoutput.AnalogOutput
            method), 17
BrokenInstall, 22
browse () (labscript_utils.qtwidgets.imageoutput.BrowseButton
            method), 19
BrowseButton (class in lab-
            script_utils.qtwidgets.imageoutput), 19
button_left () (lab-
            script_utils.qtwidgets.dragdroptab.Tab
            property), 13
button_right () (lab-
            script_utils.qtwidgets.dragdroptab.Tab
            property), 13

```

B

```

base_unit (labscript_utils.unitconversions.aom.SineAom
            check () (in module labscript_utils.memprof), 20
            attribute), 7
base_unit (labscript_utils.unitconversions.detuning.detuning
            method), 21
            attribute), 7
base_unit (labscript_utils.unitconversions.example.example1
            method), 20
            attribute), 7
base_unit (labscript_utils.unitconversions.example.example2
            method), 21
            attribute), 7
base_unit (labscript_utils.unitconversions.example.example3
            method), 21
            attribute), 7
base_unit (labscript_utils.unitconversions.linear_coil_driver.BidirectionalCoilDriver
            BiDirectionalCoilDriverSettings.Settings method), 5
            attribute), 8
base_unit (labscript_utils.unitconversions.NovaTechDDS9m.NovaTechDDS9m
            NovaTechDDS9mAmpConvex.Connection method),
            attribute), 8
base_unit (labscript_utils.unitconversions.NovaTechDDS9m.NovaTechDDS9m
            NovaTechDDS9mFreqConversion
            attribute), 8
base_unit (labscript_utils.unitconversions.optotunelens.OptotuneLens
            method), 4
            attribute), 8
base_unit (labscript_utils.unitconversions.quad_driver.quad_driver
            base_unit (labscript_utils.unitconversions.quad_monitor.quad_monitor
            attribute), 9
            attribute), 9
base_unit (labscript_utils.unitconversions.test.test_at-
            tribute), 6
before_get () (lab-
            script_utils.labconfig.EnvInterpolation
            method), 4
BG (labscript_utils.splash.Splash attribute), 11
BidirectionalCoilDriver (class in lab-
            script_utils.unitconversions.linear_coil_driver),
            8
block_combobox_signals() (lab-
            script_utils.qtwidgets.analogoutput.AnalogOutput
            method), 17
k_combobox_signals() (lab-
            script_utils.qtwidgets.enumoutput.EnumOutput
            method), 18
block_spinbox_signals() (lab-
            script_utils.qtwidgets.analogoutput.AnalogOutput
            method), 17
BrokenInstall, 22
browse () (labscript_utils.qtwidgets.imageoutput.BrowseButton
            method), 19
BrowseButton (class in lab-
            script_utils.qtwidgets.imageoutput), 19
button_left () (lab-
            script_utils.qtwidgets.dragdroptab.Tab
            property), 13
button_right () (lab-
            script_utils.qtwidgets.dragdroptab.Tab
            property), 13

```

C

```

check () (labscript_utils.filewatcher.FileWatcher
            method), 21
check () (labscript_utils.memprof.MemoryProfiler
            method), 20
check () (labscript_utils.modulewatcher.ModuleWatcher
            method), 21
check_version () (in module lab-
            script_utils.versions), 22
close () (labscript_utils.h5_lock.File method), 9
compare_to () (lab-
            script_utils.connections.ConnectionTable
            method), 5
connect_to_zclock_server () (in module lab-
            script_utils.ls_zprocess), 10

```

```

connect_value_change() (lab- attribute), 9
    script_utils.qtwidgwidgets.analogoutput.AnalogOutputderived_units (lab-
        method), 17
    script_utils.qtwidgwidgets.enumoutput.EnumOutputderived_units (lab-
        method), 18
    Connection (class in labscript_utils.connections), 3
    ConnectionTable (class in lab-
        script_utils.connections), 3
    Context (class in labscript_utils.ls_zprocess), 9
    contextMenuEvent () (lab-
        script_utils.qtwidgwidgets.imageoutput.ImageView
        method), 19
    count_types () (lab-
        script_utils.memprof.MemoryProfiler method), 20
    create_dialog () (labscript_utils.settings.Settings
        method), 5
    create_linked_width_group () (lab-
        script_utils.qtwidgwidgets.toolpalette.ToolPaletteGroup
        method), 16

D

d_MHz_from_base () (lab-
    script_utils.unitconversions.detuning.detuning
    method), 7
d_MHz_to_base () (lab-
    script_utils.unitconversions.detuning.detuning
    method), 7
data () (labscript_utils.qtwidgwidgets.dragdroptab.Tab
    property), 13
DDSSOutput (class in lab-
    script_utils.qtwidgwidgets.ddsoutput), 18
debug (class in labscript_utils.qtwidgwidgets.dragdroptab),
    14
DEBUG (labscript_utils.qtwidgwidgets.dragdroptab.debug at-
    tribute), 14
dedent () (in module labscript_utils), 3
depth (labscript_utils.qtwidgwidgets.dragdroptab.debug at-
    tribute), 14
derived_units (lab-
    script_utils.unitconversions.detuning.detuning
    attribute), 8
derived_units (lab-
    script_utils.unitconversions.linear_coil_driver.Bidirectional
    attribute), 8
derived_units (lab-
    script_utils.unitconversions.optotunelens.OptotuneLens
    attribute), 8
derived_units (lab-
    script_utils.unitconversions.quad_driver.quad_driver
    attribute), 9
derived_units (lab-
    script_utils.unitconversions.quad_monitor.quad_monitor
    attribute), 9

attribute), 9
    derived_units (lab-
        script_utils.unitconversions.test.test
        attribute), 6
    deserialise() (in module lab-
        script_utils.properties), 6
detuned_MHz_from_base () (lab-
    script_utils.unitconversions.example.example2
    method), 7
detuned_MHz_to_base () (lab-
    script_utils.unitconversions.example.example2
    method), 7
detuning (class in lab-
    script_utils.unitconversions.detuning), 7
dict_diff () (in module labscript_utils.dict_diff), 4
diff () (labscript_utils.connections.Connection
    method), 3
DigitalOutput (class in lab-
    script_utils.qtwidgwidgets.digitaloutput), 18
disable () (in module lab-
    script_utils.double_import_denier), 21
disable () (in module labscript_utils.impprof), 19
disable () (labscript_utils.impprof._ProfilingImporter
    method), 19
disconnect_value_change () (lab-
    script_utils.qtwidgwidgets.analogoutput.AnalogOutput
    method), 17
disconnect_value_change () (lab-
    script_utils.qtwidgwidgets.enumoutput.EnumOutput
    method), 18
distance_from_base () (lab-
    script_utils.unitconversions.optotunelens.OptotuneLens
    method), 8
distance_to_base () (lab-
    script_utils.unitconversions.optotunelens.OptotuneLens
    method), 8
do_update_widget_positions () (lab-
    script_utils.qtwidgwidgets.headerview_with_widgets.HorizontalHeader
    method), 15
dotdict (class in labscript_utils.testing_utils), 5
DoubleImportDenier (class in lab-
    script_utils.double_import_denier), 20
drag_in_progress () (lab-
    script_utils.qtwidgwidgets.dragdroptab.DragDropTabBar
    attribute), 12
DragDropTabBar (class in lab-
    script_utils.qtwidgwidgets.dragdroptab), 12
DragDropTabWidget (class in lab-
    script_utils.qtwidgwidgets.dragdroptab), 13
dragged_tab_index () (lab-
    script_utils.qtwidgwidgets.dragdroptab.DragDropTabBar
    property), 12
dragged_tab_parent () (lab-
    script_utils.qtwidgwidgets.dragdroptab.DragDropTabBar
    attribute), 12
dragged_tab_parent () (lab-
    script_utils.qtwidgwidgets.dragdroptab.DragDropTabBar
    attribute), 12

```

```

    property), 12
duration () (labscrip_utils.qtwidgts.dragdroptab.TabAnimation script_utils.unitconversions.example), 6
    method), 14

E
elide_label () (in module lab- example1 (class in lab-
    script_utils.qtwidgts.elide_label), 15
ElidedLabelContainer (class in lab- example2 (class in lab-
    script_utils.qtwidgts.elide_label), 14
elideMode () (labscrip_utils.qtwidgts.elide_label.ElidedLabelContainer(class in labscrip_utils.filewatcher), 21
    method), 15
filter () (labscrip_utils.setup_logging.LessThanFilter
ElideScrollArea (class in lab- method), 20
    script_utils.qtwidgts.elide_label), 14
enable () (in module lab- find_by_name () (lab-
    script_utils.double_import_denier), 21
    script_utils.connections.Connection method),
enable () (in module labscrip_utils.impprof), 19
enable () (labscrip_utils.impprof._ProfilingImporter
    method), 19
find_by_name () (lab-
ensure_connected_to_zlog () (in module lab- find_child() (lab-
    script_utils.ls_zprocess), 10
    script_utils.connections.Connection method),
ensure_running () (lab- 3
    script_utils.qtwidgts.dragdroptab.TabAnimation find_child()
    method), 14
ensure_visible () (lab- script_utils.connections.ConnectionTable
    script_utils.qtwidgts.dragdroptab.DragDropTabBar find_spec () (labscrip_utils.double_import_denier.DoubleImportDenie
    method), 12
method), 21

EnumOutput (class in lab- FingerTabBarWidget (class in lab-
    script_utils.qtwidgts.enumoutput), 18
    script_utils.qtwidgts.fingertab), 15
EnvInterpolation (class in lab- FingerTabWidget (class in lab-
    script_utils.labconfig), 4
    script_utils.qtwidgts.fingertab), 15
Event () (in module labscrip_utils.ls_zprocess), 10
FLUSH_GAP (labscrip_utils.qtwidgts.dragdroptab.DragDropTabBar
event () (labscrip_utils.elide_label.ElidedLabelContainattribute), 12
    method), 15
focusInEvent () (lab-
event () (labscrip_utils.elide_label.ElideScrollArea script_utils.qtwidgts.analogoutput.NoStealFocusDoubleSpinBox
    method), 14
method), 18
eventFilter () (lab- focusOutEvent () (lab-
    script_utils.qtwidgts.analogoutput.AnalogOutput script_utils.qtwidgts.analogoutput.NoStealFocusDoubleSpinBox
    method), 17
method), 18
eventFilter () (lab- fraction_from_base () (lab-
    script_utils.qtwidgts.digitaloutput.DigitalOutput script_utils.unitconversions.aom.SineAom
    method), 18
method), 7
eventFilter () (lab- fraction_to_base () (lab-
    script_utils.qtwidgts.enumoutput.EnumOutput script_utils.unitconversions.aom.SineAom
    method), 18
method), 7
eventFilter () (lab- Gauss_from_base () (lab-
    script_utils.qtwidgts.headerview_with_widgets.HorizontalHeaderViewWithWidgets
    method), 15
method), 6
eventFilter () (lab- Gauss_from_base () (lab-
    script_utils.qtwidgts.imageoutput.BrowseButton script_utils.unitconversions.example.example1
    method), 19
method), 7
eventFilter () (lab- Gauss_from_base () (lab-
    script_utils.qtwidgts.imageoutput.ImageOutput script_utils.unitconversions.test.test
    method), 19
method), 6
Gauss_to_base () (lab-
    script_utils.unitconversions.example.example1

```

```

        method), 7
Gauss_to_base()           (lab-      script_utils.qtwidgets.toolpalette.ToolPaletteGroup
    script_utils.unitconversions.test.test   method),  get_palette_by_index()           (lab-
    6                                         method), 16
Gcm_from_base()           (lab-      script_utils.qtwidgets.toolpalette.ToolPaletteGroup
    script_utils.unitconversions.quad_driver.quad_driver.get_sub_widget()           (lab-
    method), 9                                         method), 16
Gcm_from_base()           (lab-      script_utils.qtwidgets.ddsoutput.DDSOutput
    script_utils.unitconversions.quad_monitor.quad_monitor.convertion_class() (in module lab-
    method), 9                                         script_utils.unitconversions), 6
Gcm_to_base()             (lab-      get_value()          (labscript_utils.settings.Settings
    script_utils.unitconversions.quad_driver.quad_driver     method), 5
    method), 9                                         get_version() (in module labscript_utils.versions),
Gcm_to_base()             (lab-      23
    script_utils.unitconversions.quad_monitor.quad_monitor
    method), 9                                         H
get() (in module labscript_utils.properties), 6           h (labscript_utils.splash.Splash attribute), 11
get_AI() (labscript_utils.qtwidgets.analoginput.AnalogInput.set_locks_onto_h5py() (in module lab-
    method), 17                                         script_utils.h5_lock), 9
get_AO() (labscript_utils.qtwidgets.analogoutput.AnalogOutput.set_processor() (in module labscript_utils.ls_zprocess), 10
    method), 17                                         hardware_from_base() (lab-
getAttachedDevices()      (lab-      script_utils.unitconversions.NovaTechDDS9m.NovaTechDDS9mA
    script_utils.connections.ConnectionTable
    method), 4                                         method), 8
get_attribute() (in module lab-      hardware_to_base() (lab-
    script_utils.properties), 6                                         script_utils.unitconversions.NovaTechDDS9m.NovaTechDDS9mA
    method), 6                                         method), 8
get_attributes() (in module lab-      has_palette() (lab-
    script_utils.properties), 6                                         script_utils.qtwidgets.toolpalette.ToolPaletteGroup
    method), 6                                         method), 16
get_clean_modified_info() (lab-      hide_palette() (lab-
    script_utils.filewatcher.FileWatcher
    method), 21                                         script_utils.qtwidgets.toolpalette.ToolPaletteGroup
    method), 16                                         method), 16
get_config() (in module lab-      hide_palette_by_index() (lab-
    script_utils.ls_zprocess), 10                                         script_utils.qtwidgets.toolpalette.ToolPaletteGroup
    method), 10                                         method), 16
get_DO() (labscript_utils.qtwidgets.digitaloutput.DigitalOutput.set_sub_widget() (lab-
    method), 18                                         script_utils.qtwidgets.ddsoutput.DDSOutput
    method), 18                                         method), 16
get_EO() (labscript_utils.qtwidgets.enumoutput.EnumOutput.set_processor() (in module lab-
    method), 19                                         script_utils.qtwidgets.toolpalette.ToolPaletteGroup
    method), 19                                         method), 16
get_Image() (labscript_utils.qtwidgets.imageoutput.ImageOutput.set_processor() (in module lab-
    method), 18                                         hideSection() (lab-
    method), 19                                         script_utils.qtwidgets.toolpalette.ToolPaletteGroup
    method), 19                                         method), 15
get_import_path() (in module lab-      hideSection() (lab-
    script_utils.versions), 22                                         script_utils.qtwidgets.headerview_with_widgets.HorizontalHeader
    method), 22                                         viewWithWidgets
    method), 16                                         (class
    method), 16                                         in
get_modified_info() (lab-      HorizontalHeaderViewWithWidgets
    script_utils.filewatcher.FileWatcher
    method), 21                                         (script_utils.qtwidgets.headerview_with_widgets), 15
get_modified_times() (lab-      I_from_base() (lab-
    script_utils.filewatcher.FileWatcher
    method), 21                                         script_utils.unitconversions.optotunelens.OptotuneLens
    method), 8
get_name_from_index() (lab-      I_to_base() (labscript_utils.unitconversions.optotunelens.OptotuneLens
    script_utils.qtwidgets.toolpalette.ToolPaletteGroup
    method), 16                                         method), 8
get_palette() (lab-      icon() (labscript_utils.qtwidgets.dragdroptab.Tab
    method), 13                                         property), 13

```

```
icon_frac (labscrip_utils.splash.Splash attribute), 11
ImageOutput (class in lab-
    script_utils.qtwidgits.imageoutput), 19
imageUpdated (lab-
    script_utils.qtwidgits.imageoutput.ImageOutput
    attribute), 19
ImageView (class in lab-
    script_utils.qtwidgits.imageoutput), 19
imheight (labscrip_utils.splash.Splash attribute), 11
import_or_reload () (in module labscrip_utils), 3
imwidth (labscrip_utils.splash.Splash attribute), 12
insert_new_palette () (lab-
    script_utils.qtwidgits.toolpalette.ToolPaletteGroup
    method), 16
insertion_index_at () (lab-
    script_utils.qtwidgits.dragdroptab.DragDropTabBar
    method), 12
insertWidget () (lab-
    script_utils.qtwidgits.toolpalette.ToolPalette
    method), 16
instance () (labscrip_utils.ls_zprocess.Context class
    method), 9
instance () (labscrip_utils.ls_zprocess.ProcessTree
    class method), 10
instance () (labscrip_utils.ls_zprocess.ZMQClient
    class method), 10
InvertedDigitalOutput (class in lab-
    script_utils.qtwidgits.digitaloutput), 18
is_dragged_tab () (lab-
    script_utils.qtwidgits.dragdroptab.DragDropTabBar
    method), 12
is_json () (in module labscrip_utils.properties), 6

L
LabConfig (class in labscrip_utils.labconfig), 4
LabConfig.NoOptionError, 4
LabConfig.NoSectionError, 4
labscrip_utils
    module, 3
labscrip_utils.connections
    module, 3
labscrip_utils.dict_diff
    module, 4
labscrip_utils.double_import_denier
    module, 20
labscrip_utils.filewatcher
    module, 21
labscrip_utils.h5_lock
    module, 9
labscrip_utils.memprof
    module, 19
labscrip_utils.labconfig
    module, 4
labscrip_utils.ls_zprocess
    module, 9
labscrip_utils.modulewatcher
    module, 21
labscrip_utils.properties
    module, 6
labscrip_utils.qtwidgits
    module, 12
labscrip_utils.qtwidgits.analoginput
    module, 17
labscrip_utils.qtwidgits.analogoutput
    module, 17
labscrip_utils.qtwidgits.ddsoutput
    module, 18
labscrip_utils.qtwidgits.digitaloutput
    module, 18
labscrip_utils.qtwidgits.dragdroptab
    module, 12
labscrip_utils.qtwidgits.elide_label
    module, 14
labscrip_utils.qtwidgits.enumoutput
    module, 18
labscrip_utils.qtwidgits.fingertab
    module, 15
labscrip_utils.qtwidgits.headerview_with_widgets
    module, 15
labscrip_utils.qtwidgits.imageoutput
    module, 19
labscrip_utils.qtwidgits.inputplotwindow
    module, 16
labscrip_utils.qtwidgits.outputbox
    module, 16
labscrip_utils.qtwidgits.toolpalette
    module, 16
labscrip_utils.remote
    module, 11
labscrip_utils.settings
    module, 5
labscrip_utils.setup_logging
    module, 20
labscrip_utils.shared_drive
    module, 11
labscrip_utils.splash
    module, 11
labscrip_utils.testing_utils
    module, 5
labscrip_utils.tracelog
    module, 20
labscrip_utils.unitconversions
    module, 6
labscrip_utils.unitconversions.aom
    module, 7
labscrip_utils.unitconversions.detuning
```

```

    module, 7                               mainloop () (labscript_utils.modulewatcher.ModuleWatcher
labscript_utils.unitconversions.example      method), 21
    module, 6                               MemoryProfiler (class in labscript_utils.memprof),
labscript_utils.unitconversions.linear_coil_dri30er
    module, 8                               MHz_from_base () (lab-
labscript_utils.unitconversions.NovaTechDDS9m script_utils.unitconversions.detuning.detuning
    module, 8                               method), 7
labscript_utils.unitconversions.optotuneMHzfrom_base () (lab-
    module, 8                               script_utils.unitconversions.NovaTechDDS9m.NovaTechDDS9mF
labscript_utils.unitconversions.quad_driver   method), 8
    module, 8                               MHz_to_base () (lab-
labscript_utils.unitconversions.quad_monitor  script_utils.unitconversions.detuning.detuning
    module, 9                               method), 7
labscript_utils.unitconversions.test         MHz_to_base () (lab-
    module, 6                               script_utils.unitconversions.NovaTechDDS9m.NovaTechDDS9mF
labscript_utils.unitconversions.UnitConversionMethod), 8
    module, 7                               minimumSize () (lab-
labscript_utils.versions                   script_utils.qtwidgets.toolpalette.ToolPalette
    module, 22                             method), 16
labscript_utils.zlock                      minimumSizeHint () (lab-
    module, 11                            script_utils.qtwidgets.dragdroptab.DragDropTabBar
labscript_utils.zlog                       method), 12
    module, 11                            minimumSizeHint () (lab-
LessThanFilter      (class      in      lab- script_utils.qtwidgets.elide_label.ElidedLabelContainer
    script_utils.setup_logging), 20        method), 15
limbo (labscript_utils.qtwidgets.dragdroptab.DragDropTabBar) minimumSizeHint () (lab-
    attribute), 12
linewidhts_from_base () (lab- script_utils.qtwidgets.elide_label.ElideScrollView
    script_utils.unitconversions.detuning.detuning method), 14
linewidhts_to_base () (lab- minimumSizeHint () (lab-
    script_utils.unitconversions.detuning.detuning script_utils.qtwidgets.toolpalette.ToolPalette
    method), 8                           method), 16
    module
load () (labscript_utils.settings.Settings method), 5
load_appconfig () (in module lab- labscrip
    script_utils.labconfig), 4
Lock () (in module labscript_utils.ls_zprocess), 10
lock () (labscript_utils.qtwidgets.analogoutput.AnalogOutput method), 17
lock () (labscript_utils.qtwidgets.digitaloutput.DigitalOutput method), 18
lock () (labscript_utils.qtwidgets.enumoutput.EnumOutput method), 19
lock () (labscript_utils.qtwidgets.imageoutput.ImageOutput method), 19
log () (in module labscript_utils.tracelog), 20

M

main () (in module labscript_utils.remote), 11
main () (in module labscript_utils.zlock), 11
main () (in module labscript_utils.zlog), 11
mainloop () (labscript_utils.filewatcher.FileWatcher method), 21
    module
labscrip
    labscrip_utils, 3
    labscrip_utils.connections, 3
    labscrip_utils.dict_diff, 4
    labscrip_utils.double_import_denier,
    20
    labscrip_utils.filewatcher, 21
    labscrip_utils.h5_lock, 9
    labscrip_utils.impprof, 19
    labscrip_utils.labconfig, 4
    labscrip_utils.ls_zprocess, 9
    labscrip_utils.memprof, 20
    labscrip_utils.modulewatcher, 21
    labscrip_utils.properties, 6
    labscrip_utils.qtwidgets, 12
    labscrip_utils.qtwidgets.analoginput,
    17
    labscrip_utils.qtwidgets.analogoutput,
    17
    labscrip_utils.qtwidgets.ddsoutput,
    18
    labscrip_utils.qtwidgets.digitaloutput,
    18

```

```
labscript_utils.qtwidgets.dragdroptab,           method), 12
    12                         mousePressEvent ()                      (lab-
labscript_utils.qtwidgets.elide_label,           script_utils.qtwidgets.dragdroptab.DragDropTabBar
    14                         method), 12
labscript_utils.qtwidgets.enumoutput, mouseReleaseEvent ()          (lab-
    18
labscript_utils.qtwidgets.fingertab,             script_utils.qtwidgets.dragdroptab.DragDropTabBar
    15                         method), 12
labscript_utils.qtwidgets.headerview_with_widgets, N
    15                         NoStealFocusDoubleSpinBox (class in lab-
labscript_utils.qtwidgets.imageoutput,           script_utils.qtwidgets.analogoutput), 18
    19                         NotFound (class in labscript_utils.versions), 22
labscript_utils.qtwidgets.InputPlotWindow NovaTechDDS9mAmpConversion (class in lab-
    16
labscript_utils.qtwidgets.outputbox,             script_utils.unitconversions.NovaTechDDS9m),
    8                         NovaTechDDS9mFreqConversion (class in lab-
labscript_utils.qtwidgets.palette,               script_utils.unitconversions.NovaTechDDS9m),
    16                         8
labscript_utils.remote, 11                      NoVersionInfo (class in labscript_utils.versions), 22
labscript_utils.settings, 5
labscript_utils.setup_logging, 20
labscript_utils.shared_drive, 11
labscript_utils.splash, 11
labscript_utils.testing_utils, 5
labscript_utils.tracelog, 20
labscript_utils.unitconversions, 6
labscript_utils.unitconversions.aom, on_columnsRemoved() (lab-
    7
labscript_utils.unitconversions.detuning,       script_utils.qtwidgets.headerview_with_widgets.HorizontalHeader
    7                         method), 15
labscript_utils.unitconversions.example, on_save() (labscript_utils.settings.Settings method), 5
    6                         on_scroll_button_clicked() (lab-
labscript_utils.unitconversions.linear_coil_driver, dragdroptab.DragDropTabBar
    8                         method), 12
labscript_utils.unitconversions.NovaTechDDS9m
    8                         on_tab_moved() (lab-
labscript_utils.unitconversions.optotunelens, open_plot_window() (lab-
    8
labscript_utils.unitconversions.quad_driver,   script_utils.qtwidgets.dragdroptab.TabAnimation
    8                         method), 14
labscript_utils.unitconversions.optotunelens, OptotuneLens (class in lab-
    8                         method), 17
labscript_utils.unitconversions.quad_monitor, script_utils.unitconversions.optotunelens),
    9                         8
labscript_utils.unitconversions.test, OutputBox (class in lab-
    6
labscript_utils.unitconversions.UnitConversionBase, script_utils.qtwidgets.outputbox), 16
    7
labscript_utils.versions, 22
labscript_utils.zclock, 11
labscript_utils.zlog, 11
ModuleWatcher (class in lab- paint_tab() (labscript_utils.qtwidgets.dragdroptab.DragDropTabBar
    script_utils.modulewatcher), 21
    7                         method), 12
monkeypatch (class in labscript_utils.testing_utils), 5
mouseMoveEvent () (lab- paintEvent () (lab-
    script_utils.qtwidgets.dragdroptab.DragDropTabBar
    script_utils.qtwidgets.dragdroptab.DragDropTabBar
    15                         method), 12
script_utils.qtwidgets.fingertab.FingerTabBarWidget
    method), 15
```

```

paintEvent () (labscript_utils.splash.Splash method),
    12
path_to_agnostic () (in module lab-
    script_utils.shared_drive), 11
path_to_local () (in module lab-
    script_utils.shared_drive), 11
PlotWindow (class in lab-
    script_utils.qtwidgits.InputPlotWindow),
    16
Power_from_base () (lab-
    script_utils.unitconversions.aom.SineAom
    method), 7
Power_to_base () (lab-
    script_utils.unitconversions.aom.SineAom
    method), 7
print_details () (lab-
    script_utils.connections.Connection method),
    3
print_details () (lab-
    script_utils.connections.ConnectionTable
    method), 4
ProcessTree (class in labscript_utils.ls_zprocess), 10
profiling_import () (lab-
    script_utils.improf._ProfilingImporter
    method), 19
properties () (lab-
    script_utils.connections.Connection property),
    3

Q
quad_driver (class in lab-
    script_utils.unitconversions.quad_driver),
    8
quad_monitor (class in lab-
    script_utils.unitconversions.quad_monitor),
    9
quit_mainloop () (lab-
    script_utils.testing_utils.ThreadTestCase
    method), 5

R
register_callback () (lab-
    script_utils.settings.Settings method), 5
RemoteProcessClient () (in module lab-
    script_utils.ls_zprocess), 10
remove () (labscript_utils.qtwidgits.toolpalette.ToolPaletteGroup
    method), 16
remove_by_index () (lab-
    script_utils.qtwidgits.toolpalette.ToolPaletteGroup
    method), 17
remove_callback () (lab-
    script_utils.settings.Settings method), 5
remove_device () (lab-
    script_utils.connections.ConnectionTable
    method), 4
remove_from_linked_width_group () (lab-
    script_utils.qtwidgits.toolpalette.ToolPaletteGroup
    method), 17
reorder_palette () (lab-
    script_utils.qtwidgits.toolpalette.ToolPaletteGroup
    method), 17
reorder_palette_by_index () (lab-
    script_utils.qtwidgits.toolpalette.ToolPaletteGroup
    method), 17
resizeEvent () (lab-
    script_utils.qtwidgits.toolpalette.ToolPalette
    method), 16
run () (labscript_utils.qtwidgits.InputPlotWindow.PlotWindow
    method), 16
run () (labscript_utils.testing_utils.ThreadTestCase
    method), 5
run_script_as_main () (lab-
    script_utils.testing_utils.ThreadTestCase
    method), 5

S
save_appconfig () (in module lab-
    script_utils.labconfig), 4
SCROLL_BUTTON_GAP (lab-
    script_utils.qtwidgits.dragdroptab.DragDropTabBar
    attribute), 12
SCROLL_BUTTON_WIDTH (lab-
    script_utils.qtwidgits.dragdroptab.DragDropTabBar
    attribute), 12
sectionSizeFromContents () (lab-
    script_utils.qtwidgits.headerview_with_widgets.HorizontalHeader
    method), 15
selected_index () (lab-
    script_utils.qtwidgits.enumoutput.EnumOutput
    property), 19
selected_option () (lab-
    script_utils.qtwidgits.enumoutput.EnumOutput
    property), 19
selected_unit () (lab-
    script_utils.qtwidgits.analogoutput.AnalogOutput
    property), 17
serialise () (in module labscript_utils.properties), 6
set_AI () (labscript_utils.qtwidgits.analoginput.AnalogInput
    method), 17
set_AO () (labscript_utils.qtwidgits.analogoutput.AnalogOutput
    method), 17
set_attributes () (in module lab-
    script_utils.properties), 6
set_combobox_model () (lab-
    script_utils.qtwidgits.analogoutput.AnalogOutput
    method), 17
set_combobox_model () (lab-
    script_utils.qtwidgits.enumoutput.EnumOutput
    method), 17

```

```

        method), 19
set_device_properties() (in module lab- script_utils.qtwidg script_utils.qtwidg script_utils.qtwidg
    script_utils.properties), 6   script_widgets.dragdroptab.DragDropTabWidget
                                setWidget () (labscript_utils.qtwidg script_utils.elide_label.ElideScrollArea
set_DO() (labscript_utils.qtwidg script_utils.digitaloutput.DigitalOutput   method), 14
                                setWidget () (labscript_utils.qtwidg script_utils.headerview_with_widgets.Horiz
    method), 18                                method), 15
set_EO() (labscript_utils.qtwidg script_utils.enumoutput.EnumOutput   method), 15
                                show () (labscript_utils.splash.Splash method), 12
set_Image() (labscript_utils.qtwidg script_utils.imageoutput.ImageOutput   palette () (lab-
    method), 19                                script_utils.qtwidg script_utils.toolpalettes.ToolPaletteGroup
set_limits() (lab-                                method), 17
    script_utils.qtwidg script_utils.analogoutput.AnalogOutput show_palette () (lab-
    method), 17                                script_utils.qtwidg script_utils.toolpalettes.ToolPaletteGroup
set_num_decimals() (lab-                                method), 17
    script_utils.qtwidg script_utils.analogoutput.AnalogOutput show_sub_widget () (lab-
    method), 17                                script_utils.qtwidg script_utils.ddsoutput.DDSOutput
set_selected_unit() (lab-                                method), 18
    script_utils.qtwidg script_utils.analogoutput.AnalogOutput showEvent () (labscript_utils.qtwidg script_utils.headerview_with_widgets.Horiz
    method), 17                                method), 15
set_spinbox_value() (lab-                                showSection () (lab-
    script_utils.qtwidg script_utils.analogoutput.AnalogOutput                                script_utils.qtwidg script_utils.headerview_with_widgets.HorizontalHeade
    method), 18                                method), 15
set_step_size() (lab- SineAom (class in labscript_utils.unitconversions.aom),
    script_utils.qtwidg script_utils.analogoutput.AnalogOutput                                7
    method), 18                                sizeHint () (labscript_utils.qtwidg script_utils.dragdroptab.DragDropTabBar
set_tab_parent() (lab-                                method), 13
    script_utils.qtwidg script_utils.dragdroptab.DragDropTabBar sizeHint () (labscript_utils.qtwidg script_utils.elide_label.ElidedLabelContainer
    method), 12                                method), 15
set_value() (labscript_utils.qtwidg AnalogInput   (labscript_utils.qtwidg script_utils.elide_label.ElideScrollArea
    method), 17                                method), 14
setElideMode() (lab-                                sizeHint () (labscript_utils.qtwidg script_utils.toolpalettes.ToolPalette
    script_utils.qtwidg script_utils.dragdroptab.DragDropTabWidget   method), 16
    method), 13                                socket () (labscript_utils.ls_zprocess.Context
setElideMode() (lab-                                method), 9
    script_utils.qtwidg script_utils.elide_label.ElidedLabelContainer start () (in module labscript_utils.memprof), 20
    method), 15                                start () (labscript_utils.memprof.MemoryProfiler
setElideMode() (lab-                                method), 14
    script_utils.qtwidg script_utils.elide_label.ElideScrollArea state () (labscript_utils.qtwidg script_utils.digitaloutput.DigitalOutput
    method), 14                                property), 18
setSectionHidden() (lab-                                property), 18
    script_utils.qtwidg script_utils.headerview_with_widgets.HorizontalHeaderWithWidgets digitaloutput.InvertedDigitalOutput
    method), 15                                property), 18
setStyleSheet() (lab- stop () (labscript_utils.filewatcher.FileWatcher
    script_utils.qtwidg script_utils.headerview_with_widgets.HorizontalHeaderWithWidgets
    method), 15                                stylesheets (labscript_utils.qtwidg script_utils.headerview_with_widgets.HorizontalHeaderWithWidgets
                                                attribute), 16
Settings (class in labscript_utils.settings), 5 (lab-
setTopLevelWindow() (lab- T
    script_utils.qtwidg InputPlotWindow.PlotWindow Tab (class in labscript_utils.qtwidg script_utils.dragdroptab), 13
    method), 16                                tab_widgets (labscript_utils.qtwidg script_utils.dragdroptab.DragDropTabBar
set_logging() (in module lab-                                attribute), 13
    script_utils.setup_logging), 20
setUsesScrollButtons() (lab- TabAnimation (class in lab-
    script_utils.qtwidg script_utils.dragdroptab.DragDropTabBar script_utils.qtwidg script_utils.dragdroptab), 14
    method), 12                                tabAt () (labscript_utils.qtwidg script_utils.dragdroptab.DragDropTabBar
setUsesScrollButtons() (lab-                                method), 13
                                                method), 13

```

```

tabInserted()                               (lab- UnidirectionalCoilDriver (class in lab-
    script_utils.qtwidgets.dragdroptab.DragDropTabBar      script_utils.unitconversions.linear_coil_driver),
    method), 13                                8

tabInserted()                               (lab- unit_conversion_params()          (lab-
    script_utils.qtwidgets.dragdroptab.TabAnimation       script_utils.connections.Connection property),
    method), 14                                3

tabLayoutChange()                          (lab- unit_list (labscript_utils.unitconversions.UnitConversionBase.UnitCon-
    script_utils.qtwidgets.dragdroptab.DragDropTabBar     attribute), 7
    method), 13

tabRect () (labscript_utils.qtwidgets.dragdroptab.DragDropTabBar script_utils.unitconversions.UnitConversionBase),
    method), 13                                7

tabRemoved()                               (lab- unload () (labscript_utils.modulewatcher.ModuleWatcher
    script_utils.qtwidgets.dragdroptab.DragDropTabBar     method), 21
    method), 13                                  unlock () (labscript_utils.qtwidgets.analogoutput.AnalogOutput

tabRemoved()                               (lab- method), 18
    script_utils.qtwidgets.dragdroptab.TabAnimation unlock () (labscript_utils.qtwidgets.digitaloutput.DigitalOutput
    method), 14                                method), 18

tabSizeHint()                            (lab- unlock () (labscript_utils.qtwidgets.enumoutput.EnumOutput
    script_utils.qtwidgets.fingertab.FingerTabBarWidget   method), 19
    method), 15                                unlock () (labscript_utils.qtwidgets.imageoutput.ImageOutput

target () (labscript_utils.qtwidgets.dragdroptab.TabAnimation method), 19
    method), 14                                update () (labscript_utils.qtwidgets.dragdroptab.DragDropTabBar
                                                method), 13

tau (labscript_utils.qtwidgets.dragdroptab.TabAnimation
    attribute), 14                           update_dragged_tab_animation_pos () (lab-
                                                script_utils.qtwidgets.dragdroptab.DragDropTabBar
    method), 13

test (class in labscript_utils.unitconversions.test), 6
text () (labscript_utils.qtwidgets.dragdroptab.Tab
    property), 13                           update_elide_widget ()          (lab-
                                                script_utils.qtwidgets.elide_label.ElidedLabelContainer
    method), 15

thinspace (labscript_utils.qtwidgets.headerview_with_widgets.HorizontalHeaderViewWithWidgetWatcher
    attribute), 16                           update_files ()           (lab-
                                                script_utils.qtwidgets.dragdroptab.DragDropTabBar
    method), 21

ThreadTestCase (class in lab-
    script_utils.testing_utils), 5           update_indent ()          (lab-
                                                script_utils.qtwidgets.headerview_with_widgets.HorizontalHeader-
    method), 16

ToolPalette (class in lab-
    script_utils.qtwidgets.toolpalette), 16
ToolPaletteGroup (class in lab-
    script_utils.qtwidgets.toolpalette), 16
tooltip () (labscript_utils.qtwidgets.dragdroptab.Tab
    property), 14                           update_scroll_button_state () (lab-
                                                script_utils.qtwidgets.dragdroptab.DragDropTabBar
    method), 13

trace () (labscript_utils.qtwidgets.dragdroptab.debug
    class method), 14                         update_tab_index ()          (lab-
                                                script_utils.qtwidgets.dragdroptab.DragDropTabBar
    method), 13

U

unblock_combobox_signals () (lab-
    script_utils.qtwidgets.analogoutput.AnalogOutput
    method), 18                           update_text ()          (labscript_utils.splash.Splash
                                                method), 12

unblock_combobox_signals () (lab-
    script_utils.qtwidgets.enumoutput.EnumOutput
    method), 19                           update_widget_positions () (lab-
                                                script_utils.qtwidgets.headerview_with_widgets.HorizontalHeader-
    method), 16

unblock_spinbox_signals () (lab-
    script_utils.qtwidgets.analogoutput.AnalogOutput
    method), 18                           updateCurrentTime () (lab-
                                                script_utils.qtwidgets.dragdroptab.TabAnimation
    method), 14

updateMinimumSize () (lab-
    script_utils.qtwidgets.toolpalette.ToolPalette

```

method), 16
zmq_push_string() (in module labscript_utils.ls_zprocess), 10

V
value () (labscript_utils.qtwidgits.imageoutput.ImageOutput property), 19
vectorise() (in module labscript_utils.unitconversions.UnitConversionBase), 7

VersionException, 22
viewportEvent() (labscript_utils.qtwidgits.headerview_with_widgets.HorizontalHeaderViewWithWidgets method), 16

W

w (labscript_utils.splash.Splash attribute), 12
W_from_base() (labscript_utils.unitconversions.example.example3 method), 7

W_to_base() (labscript_utils.unitconversions.example.example3 method), 7

wait_for() (labscript_utils.testing_utils.ThreadTestCase static method), 5

whats_this() (labscript_utils.qtwidgits.dragdroptab.Tab property), 14

wheelEvent() (labscript_utils.qtwidgits.analogoutput.NoStealFocusDoubleSpinBox method), 18

widget() (labscript_utils.qtwidgits.dragdroptab.Tab property), 14

widgetAt() (labscript_utils.qtwidgits.dragdroptab.DragDropTabBar method), 13

widths_linked() (labscript_utils.qtwidgits.toolpalette.ToolPaletteGroup property), 17

write_to_file() (labscript_utils.memprof.MemoryProfiler method), 20

Z

zmq_get() (in module labscript_utils.ls_zprocess), 10
zmq_get_multipart() (in module labscript_utils.ls_zprocess), 10

zmq_get_raw() (in module labscript_utils.ls_zprocess), 10

zmq_get_string() (in module labscript_utils.ls_zprocess), 10

zmq_push() (in module labscript_utils.ls_zprocess), 10

zmq_push_multipart() (in module labscript_utils.ls_zprocess), 10

zmq_push_raw() (in module labscript_utils.ls_zprocess), 10